

# Kouman pou w sèvi ak inalatè a avèk espasè a ak mask la



1. Retire bouchon aspiratè a



2. Souke inalatè a



3. Tache l ak espasè a



4. Chita dwat; kouvri nen ak bouch ou avèk mask la



5. Peze isit la. Sa fè medikaman an desann nan espasè a



6. Rale souf ou antre ak pouse l sòti sis fwa



Si w bezwen yon lòt boufe medikaman an, tann 1 minit epi refè etap 4 jiska 6 yo



7. Rense bouch ou ak dlo epi **krache**

