

# How to use your nebulizer



Neighborhood  
Health Plan™



1. Wash your hands



2. Put nebulizer machine on flat, steady surface



3. Plug into outlet



4. Connect one end of tubing to the nebulizer



5. Add medication



6. Attach mouthpiece to the medication holder



7. Connect the other end of the tubing to the medication holder



8. Turn on the nebulizer machine



WRONG! Do NOT hold mouthpiece in front of face



9. Put lips around mouthpiece. Take **slow, deep** breaths through your **MOUTH**



10. When finished, wash parts in soapy water and rinse



11. Allow parts to air dry