

Kouman pou w sèvi ak inalatè a avèk OptiHaler®




1. Retire kanistè la nan kontènè plastik la



2. Retire kouvèti pyès pou bouch la



3. Mete kanistè la nan OptiHaler la konsa



4. Veye byen pou kouvèti nan bout la fèmen byen sere



5. Souke byen



6. Pouse tout souf ou sòti



7. Kole bouch ou nan pyès pou bouch la



8. Pese desann **AN MENM TAN** w ap rale souf ou **ANTRE** dousman jouk nan fon



9. Kenbe souf ou pandan 10 segonn si w kapab



10. Pouse souf ou sòti dousman



Si w bezwen yon lòt boufe medikaman an...



...tann **YON** minit, epi refè etap 8 pou rive 10 yo



11. Remete bouchon an



12. Rense bouch ou ak dlo epi **krache** dlo'a

